Flocculator

SW6


Instructions for use

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
Fig 1 SW6


Fig 2 Control panel

Thank you for purchasing this Bibby Scientific product. To get the best performance from the equipment, and for your personal safety, please read these instructions carefully before use.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuations not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

Electrical Installation

THIS EQUIPMENT MUST BE EARTHED

Before connection please ensure that the line supply corresponds to that stated on the rating label.

Power requirements:

Model	Voltage & Wattage
SW6	230V, 50Hz, ~ , 200W

There is an IEC socket at the side of the instrument for connection to the mains supply. The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. See the enclosed instruction sheet for advice on how to carry out this procedure.

Should the mains lead require replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC 320 plug should be used.

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

The mains lead should be connected to the instrument BEFORE connection to the mains supply.

Safety Advice Before Use

- ❖ Never move or carry the unit when in use or connected to the mains electricity supply.
- ❖ Do not mix flammable liquids or use the equipment in hazardous atmospheres.
- ❖ In the case of mains interruption, a fault or mechanical failure, the unit will NOT continue to operate on removal of fault.
- ❖ Mechanical energy can lead to breakage of glass vessels. Use with care.

General Description

The Stuart SW6 Flocculator is designed to ensure repeatable conditions between up to six samples during water and sewage flocculation tests.

The speed control has electronic feedback to accurately maintain the rotational speed of the stainless steel paddles. Digital displays clearly indicate the speed, which is variable between 25 and 250 rpm and time, which can be set to count down from 1 to 99 minutes. Alternatively, the use of two adjustable pre-set programs (P1 and P2) allows quick and accurate repetition of pre-selected speeds and times.

All paddles can be easily lifted out of the samples when required. Observation and monitoring of samples is assisted by diffused base illumination and a black background.

One litre Pyrex® glass beakers and a coagulant dispenser are available as accessories.

Operation (see Fig 1 & 2)

Place the unit on a firm, level surface and connect to the mains supply. Leave sufficient space to access the side IEC socket.

Turn the unit ON at the side mains On / Off switch and observe the fluorescent lamp illuminate.

To position a beaker, lift the paddle up to its maximum height using the paddle lift and place the vessel over the illuminated base. Once all the samples are in position, lower the paddles into the beakers to the required depth.

Selecting a speed / time

Press the <MAN> button to selected the manual mode. To set the speed, hold down the <SET> key under the RPM display and use the <UP> and <DOWN> arrow keys until the desired speed appears on the digital display.

To set a running time, hold down the <SET> key under the TIME display and use the <UP> and <DOWN> arrow keys until the desired running time in minutes appears on the digital display.

IF A RUNNING TIME IS NOT REQUIRED, SET 00 IN THE TIME DISPLAY.

To start the run press the <START / STOP> key. The paddles will turn at the selected speed. If a time has been set, the timer will count down.

When time zero is reached, the unit will beep and the paddles stop rotating. The timer display will then flash and count up, in minutes, to indicate time elapsed since the paddles stopped.

The program can be halted at any time by pressing the <START / STOP> key. To change the settings, stop rotation and repeat the above procedure.

Pre-set programs

In order to input a pre-set speed and time, press and hold the <P1> key for 5 seconds. The green P1 light will flash. To set the speed, hold down the <SET> key under the RPM display and use the <UP> and <DOWN> arrow keys until the desired speed appears on the digital display. To set a running time, hold down the <SET> key under the TIME display and use the <UP> and <DOWN> arrow keys until the desired time in minutes appears on the digital display.

IF A RUNNING TIME IS NOT REQUIRED, SET 00 IN THE TIME DISPLAY.

Once the correct speed and time have been set, press the <P1> key to confirm. The settings will be held in the memory.

Pressing the <P1> key at any time during operation will select the preset values. A green LED shows that the P1 program has been selected.

If another preset speed and time is required press and hold down the <P2> key and repeat the above setting procedure.

New P1 and P2 values can be inputted at any time by repeating the above procedure.

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

The following spares and accessories are available from your laboratory dealer.

Description	Catalogue Number
Pyrex® glass beakers, (pk. 10)	1015/22D
Coagulant Dispenser	SW6/2

For a comprehensive list of parts required by service engineers conducting internal repairs and a service manual, please contact the Service Department at Bibby Scientific Ltd. quoting both the model and serial number.

Only spare parts supplied or specified by Bibby Scientific Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, please contact the Service Department of Bibby Scientific Ltd. or the point of sale.

Warranty

Bibby Scientific Ltd warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim, Bibby Scientific Ltd will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorized personnel.

Bibby Scientific Ltd.
Beacon Road
Stone, Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0) 1785 812121
Fax: +44 (0) 1785 813748
e-mail info@bibby-scientific.com
www.bibby-scientific.com

Technical Specification

	SW6
No. of samples	6
Speed range	25-250rpm
Timer	0 to 99 minutes
Pre-set programs	2
Digital displays	LED
Dimensions	750 x 210 x 460mm
Net weight	17kg
Power	200W


Fig 1 SW6


Fig 2 Tableau de commande

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Notes

These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards and we cannot guarantee that interference will not occur

in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Service Department of Bibby Scientific Ltd.

Declaration of Conformity

Model SW6 FLOCCULATOR

This product complies with the requirements of the EU Directives listed below:

89/336/EEC Electromagnetic Compatibility Directive, amended by 92/31/EEC and 93/68/EEC.

73/23/EEC Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN61326: 1997, EN61000-4-2: 1995, EN61000-4-4: 1995, EN61000-4-5: 1995
EN61000-4-6: 1996, EN61000-4-11: 1994, EN61000-3-2: 1995, EN61000-3-3: 1995
EN61010-1: 2001, EN61010-2-051: 1996.

The technical documentation required to demonstrate that the products meet the requirements of these Standards is available for inspection on request.

Signed:  (Mr C. R. Baggaley)

Date: 22/8/01

Authority: Technical Manager

for

Bibby Sterilin Ltd - Stone - Staffs - ST15 OSA - UK
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

INSPECTION REPORT

MODEL SW6

ELECTRICAL SAFETY

1. Earth continuity
2. Insulation
3. Flash test


FUNCTIONAL

1. Indicators
2. Speed / time
3. Visual acceptance

QUALITY CONTROL INSPECTOR

Bibby Scientific France SAS

ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bsf@bibby-scientific.fr
www.bibby-scientific.com

Bibby Scientific Italia Srl

Via Alcide de Gasperi 56
20077 Riozzo di Cerro al Lambro
Milano Italia
Tel: +39 (0)2 98230679
Fax: +39 (0)2 98230211
e-mail: marketing@bibby-scientific.it
www.bibby-scientific.it

Bibby Scientific Ltd

Beacon Road Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
e-mail: info@bibby-scientific.com
www.bibby-scientific.com

Bibby Scientific US Ltd

3 Terri Lane Suite 10
Burlington NJ 08016
USA
Tel: 800-225-9243
Fax: 609-589-2571
www.bibby-scientific.com

Bibby Scientific (Asia) Ltd

Room 607 Yen Sheng Centre
64 Hoi Yuen Road Kwun Tong
Kowloon Hong Kong
Tel: +852 3583 1581
Fax: +852 3583 1580
e-mail bibby@bibby-scientificasia.com
www.bibby-scientificasia.com