



A

DURAN® YOUTILITY Glass Bottles



Complete DURAN YOUTILITY bottle sets. GL45 screw caps, pouring rings, labels and bottle tags.

- Duran borosilicate glass 3.3
- Ergonomically shaped bottle and screw cap
- Colourful bottle tags
- Self adhesive, bottle and cap labels
- Performance range: -40°C to +150°C
- Includes autoclave/steam sterilisation Class 1 process indicator (ANSI/AAMI/ISO 11140-1:2005)

Supplied in packs of 4

Code	Alt Ref	Description	Capacity, mL
BOT1029	1502527	Clear glass	125
BOT1030	1502528	Clear glass	250
BOT1031	1502534	Clear glass	500
BOT1032	1502535	Clear glass	1000
BOT1033	1525322	Amber glass	125
BOT1034	1525324	Amber glass	250
BOT1035	1525325	Amber glass	500
BOT1036	1525356	Amber glass	1000

Accessories

Code	Alt Ref	Description	Pack
BOT1037	1503039	Replacement GL45 screw caps	10
BOT1038	1503041	Replacement pouring rings	16
BOT1039	1503043	GL45 bottle tags, mixed colours	16
BOT1041	1525320	Self-adhesive label sets, white polyester	200



NEW

DURAN TILT Glass Media Bottles



DURAN TILT - a smart design for a cell media preparation bottle. Designed for ease of handling and pipetting in biosafety cabinets and under clean hoods:

- Non-cytotoxic DURAN glass
- Wide neck for serological pipettes
- Vacuum resistant
- Upright and tilted positions - dual scale

Code	Alt Ref	Description	Pack
BOT1100	218914459	Bottle clear glass with GL56 polypropylene screw cap 500mL	4

Accessories

Code	Alt Ref	Description	Pack
BOT1102	29 243 56 01	Light shield and 4x GL56 bottle tags (orange - yellow - blue - purple)	4
BOT1104	29 243 56 26	GL56 bottle tag orange silicone	20
BOT1106	29 243 56 34	GL56 bottle tag yellow silicone	20
BOT1108	29 243 56 59	GL56 bottle tag blue silicone	20
BOT1110	29 243 56 67	GL56 bottle tag purple silicone	20
BOT1112	29 119 56 01	Thread adapter - GL56 bottle to GL45 fitting	Each
BOT1114	29 401 56 04	Labels for GL56 cap	60
BOT1116	29 229 56 02	Screw cap GL56 white polypropylene	10



NEW