

## DURAN Wide Neck Media Bottles

DURAN wide neck media bottles with cap and pouring ring, made from DURAN borosilicate glass 3.3 to ISO 3585:1998. Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service; very resistant to chemical attack. Retrace code on each bottle allows downloading of individual batch certificate.

- Wide neck for scooping, pouring and stirring
- Ideal for liquids, pastes and powders
- Supplied with quick release GLS80 PP cap and pouring ring



Code	Alt Ref	Capacity, mL	Pack
<b>BOT2117</b>	218603656	250	10
<b>BOT2110</b>	1112627	500	10
<b>BOT2112</b>	1112713	1000	10
<b>BOT2114</b>	1112715	2000	10
<b>BOT2111</b>	218606953	3500	Each
<b>BOT2116</b>	1113949	5000	Each
<b>BOT2118</b>	1113950	10,000	Each
<b>BOT2120</b>	1113951	20,000	Each
<b>BOT2124</b>	1200265	30,000	Each
<b>BOT2126</b>	1200154	50,000	Each
<b>BOT2122</b>	1112716	Spare cap GLS80	10
<b>BOT2119</b>	1160166	Spare pouring ring GLS80	10



## DURAN® Baffled Wide Mouth Bottles GLS 80®

The DURAN® GLS 80® baffled, wide mouth bottle has three vertical baffles into the inner surface, the swirling motion is disrupted, and a much more effective mixing occurs. The baffles improve the top-to-bottom circulation and produce a greater radial mixing. The baffles also offer other advantages such as suppressing the formation of a central vortex, and reducing dead-zones where mixing does not occur.

- Compatible with DURAN® GLS 80® magnetic stirrers, other magnetic stir bars, and with the full range of wide mouth GLS 80® connection cap systems for liquid addition or removal when mixing with stir bars
- The baffles are moulded into the glass side wall of the bottles and form a continuous, smooth inner surface making clean up easy
- Bottles are fully autoclavable
- Easy-to-read graduated volume scale
- Manufactured from DURAN® borosilicate 3.3 glass

Code	Alt Ref	Description
<b>BOT1132</b>	212863658	250mL
<b>BOT1134</b>	212864457	500mL
<b>BOT1136</b>	212865453	1000mL

## Stirred Reactors

DURAN® Stirred reactor wide mouth cap GLS 80® supplied with magnetic stirrer, 4 connection ports (2x GL18 and 2x GL14).

**Note:** Does not include red PBT caps as shown in image.

Code	Alt Ref	Description
<b>BOT1138</b>	212639107	Suitable for use with 250mL Bottles
<b>BOT1140</b>	1200379	Suitable for use with 500mL or 1000mL Bottles