



Output Up to 80L/hr
Purity Type 3
Feedwater Pretreated

NEW

SLS Lab Pro PurA-Q Primo Type 3 Water Purification Systems



Using the most advanced Reverse Osmosis (RO) technology available, the PURA-Q Primo efficiently removes up to 99% of the contaminants from your mains water feed to produce Primary Grade Type 3 RO water. With a choice of four flow rates up to 80L/hr and a variety of storage options, feeding the critical equipment in your laboratory with RO water has never been so made-to-measure.

- Small footprint solution to feed glass washers or other laboratory equipment
- >50% recovery of RO water
- Low running costs
- Excellent monitoring of conductivity
- One easy to change pre-filter
- Make-up rates ranging from 10-80 litres per hour
- Bench, tank or wall mounting options
- Fed directly from your mains water supply

Specifications	
Inorganics Rejection Rate, %	98
Feed Water Pressure, bar	0.1 to 6 bar
TOC, ppb	≤5
Dims, w x d x h, mm	490 x 290 x 500
Weight, kg	23

Code	Model	Make-Up Rate, L/hr	Price
WAT2006	Primo 10	10	£4163.00
WAT2008	Primo 20	20	£4541.00
WAT2016	Primo 50	30	£6433.00
WAT2018	Primo 80	40	£7194.00

Reservoirs

Code	Capacity, L	Description	Dims, w x d x h, mm	Price
WAT2444	30	Storage reservoir	405 x 230 x 495	£780.00
WAT2446	60	Storage reservoir	480 x 330 x 495	£913.00
WAT2168	100	Storage reservoir	480 x 360 x 810	£1340.00
WAT2190	100	Storage reservoir with pump	480 x 360 x 810	£2255.00

Accessories and Consumables

Code	Description	Price
WAT2178	RO membrane for 20L/hr	£567.00
WAT2180	RO membrane for 10L/hr	£509.00
WAT2182	RO membrane for 50 & 80L/hr	£624.00
WAT2172	Valve kit for washing machine feed (single POU)	£151.00
WAT2452	Wall mounting kit	£112.00
WAT2448	Wall mounting kit for 30L tank	£66.00
WAT2450	Wall mounting kit for 60L tank	£80.00
WAT2174	UV upgrade kit for 30 & 60L tank	£653.00
WAT2176	UV upgrade kit for 100L tank	£806.00
WAT2423	Pre-treatment module	£168.00
WAT2430	Air vent filter for any tank	£170.00