

## Sartorius arium® Systems Cont.



### Sartorius arium® pro

Offers a flexible and modular system which is cost effective. All systems meet and exceed the ASTM Type 1 water quality standards and can dispense up to 2L/min of consistently high-quality ultrapure water with a conductivity of 0.055µS/cm (18.2Meg-ohm/cm) each hour which is virtually microorganism-free when a Sartopore® 2 150 end filter is used. The patented Sartorius technology, SD cardslot, long service life, and low maintenance requirement mean that the arium® pro systems are easy-to-use, efficient and reliable whilst the modular approach means that systems can be tailored to specific applications.

- Options of UV and Ultrafiltration
- Cartridge packs available for analytical, biological, standard and feedwater applications
- Option of a TOC monitor



WAT4242

Output	Purity Type	Feedwater
up to 2L/min	I	Pretreated

1 YEAR



WAT4260

Output	Purity Type	Feedwater
up to 2L/min	I	Pretreated

1 YEAR

Code	Alt Ref	Description
WAT4242	H2obasicT	Basic with reduced operating features Type 1 ion exchange water system bench mounted
WAT4243	H2obasicB	Basic with reduced operating features Type 1 ion exchange water system wall mounted
WAT4244	H2oproDIT	DI Type 1 ion exchange, bench mounted
WAT4245	H2oproDIB	DI Type 1 ion exchange, wall mounted
WAT4246	H2oproDID	DI Type 1 ion exchange, under bench mounted
WAT4247	H2oproUFT	UF Type 1 ion exchange and UF, bench mounted
WAT4248	H2oproUFB	UF Type 1 ion exchange and UF, wall mounted
WAT4249	H2oproUFD	UF Type 1 ion exchange and UF, under bench mounted
WAT4250	H2oproUVT	UV Type 1 ion exchange and UV, bench mounted
WAT4251	H2oproUVB	UV Type 1 ion exchange and UV, wall mounted
WAT4252	H2oproUVD	UV Type 1 ion exchange and UV, under bench mounted
WAT4253	H2oproUVTTOC	UV Type 1 ion exchange and UV, with TOC monitor bench mounted
WAT4254	H2oproUVBTOC	UV Type 1 ion exchange and UV, with TOC monitor wall mounted
WAT4255	H2oproUVDTOC	UV Type 1 ion exchange and UV, with TOC monitor under bench mounted
WAT4256	H2oproVFT	VF Type 1 ion exchange and UV and UF, bench mounted
WAT4257	H2oproVFB	VF Type 1 ion exchange and UV and UF, wall mounted
WAT4258	H2oproVFD	VF Type 1 ion exchange and UV and UF, under bench mounted
WAT4259	H2oproVFTTOC	VF Type 1 ion exchange and UV and UF, with TOC monitor bench mounted
WAT4260	H2oproVFBTOC	VF Type 1 ion exchange and UV and UF, with TOC monitor wall mounted
WAT4261	H2oproVFDTOC	VF Type 1 ion exchange and UV and UF, with TOC monitor under bench mounted