

Thermo Scientific™ Nalgene™ Wide Mouth Bottles, PPCO

Thermo
SCIENTIFIC

- Store specimens where autoclaving is required
- Lab quality polypropylene copolymer bottles
- Excellent chemical resistance to most acids, bases and alcohol
- Translucent for easy viewing of liquid levels
- Withstands temperatures up to +121°C
- Linerless caps for leakproof performance

Code	Alt Ref	Capacity, mL	Closure Size, mm	Dims, Ø x h, mm	Pack	Price
BOT0100	2105-0001	30	28	34 x 63	12	£12.51
BOT0102	2105-0002	60	28	39 x 86	12	£14.56
BOT0104	2105-0004	125	38	50 x 99	12	£22.70
BOT0106	2105-0008	250	43	61 x 131	12	£29.70
BOT0108	2105-0016	500	53	73 x 168	12	£50.00
BOT0110	2105-0032	1000	63	91 x 199	6	£29.50



Thermo Scientific™ Nalgene™ Wide Mouth Bottles, LDPE

Thermo
SCIENTIFIC

- Translucent low density polyethylene
- Ideal for storing chemicals and standards, as well as collecting and storing distilled water
- Extremely low trace metal content in LDPE makes it ideal for trace metal analysis
- Excellent chemical resistance to most acids, bases and alcohols
- Wide mouth is easy to fill
- Linerless caps ensure leakproof performance without a liner that can wrinkle, cause leaks and contaminate solutions
- Can be used in freezer to -100°C [-148°F] for long term storage
- Also available sterile and in bulk packs - please enquire for details

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT3542	2103-0002	60	39 x 86	28	12	£13.62
BOT3544	2103-0004	125	50 x 99	38	12	£18.89
BOT3546	2103-0008	250	61 x 131	43	12	£24.10
BOT3548	2103-0016	500	73 x 168	53	12	£39.00
BOT3550	2103-0032	1000	91 x 199	63	6	£26.20



Thermo Scientific™ Nalgene™ Large Wide Mouth Bottles, HDPE

Thermo
SCIENTIFIC

- Excellent chemical resistance to most acids, bases and alcohols make this a versatile bottle
- Can be used in freezer to -100°C [-148°F] for long term storage
- Linerless cap closes the bottle without a liner that can wrinkle, cause leaks and contaminate your solutions
- Not liquid leak proof
- 100mm polypropylene screw caps

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Closure Size, mm	Price
BOT0120	2120-0005	2	119 x 244	100	£15.34
BOT0122	2120-0010	4	153 x 288	100	£21.50

