

Thermo Scientific™ Nalgene™ Narrow Mouth Bottles, HDPE

Thermo
SCIENTIFIC

- Semi-rigid, translucent all-purpose bottles
- Specifically designed for long-term, demanding lab use
- Offer excellent chemical resistance
- Versatile and indispensable in the lab
- Uniform walls provide durability and resistance to splitting or punctures
- Can be used in a freezer up to -100°C (-149°F)
- Linerless polypropylene closures provide the ultimate in leakproof protection without the use of liner that can wrinkle, cause leaks or contaminate your reagents
- Leak proof



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT0040	2002-9125	4	16 x 41	13	12	£8.57
BOT0042	2002-9025	8	25 x 45	20	12	£8.54
BOT0044	2002-9050	15	25 x 58	20	12	£9.16
BOT0046	2002-0001	30	34 x 61	20	12	£10.22
BOT0048	2002-0002	60	85 x 37	20	12	£12.30
BOT0050	2002-0004	125	50 x 101	24	12	£14.49
BOT0054	2002-0008	250	61 x 133	24	12	£21.30
BOT0056	2002-0016	500	73 x 170	28	12	£30.70
BOT0060	2002-0032	1000	91 x 216	38	6	£25.20

Thermo Scientific™ Nalgene™ Narrow Mouth Bottles, LDPE

Thermo
SCIENTIFIC

- Extremely low in trace metal content
- Excellent chemical resistance to most acids, bases and alcohols
- Guaranteed leak proof
- Can be used in freezer to -100°C (-148°F) for long-term storage
- Translucence allows easy viewing of liquid levels
- Linerless caps provide the ultimate in leakproof protection without use of liner that can wrinkle, cause leaks or contaminate your reagents
- Also available sterile and in bulk packs - please enquire for details



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT3220	2003-9025	8	25 x 45	20	12	£8.54
BOT3222	2003-9050	15	25 x 58	20	12	£9.16
BOT3224	2003-0001	30	34 x 61	20	12	£10.22
BOT3226	2003-0002	60	37 x 85	20	12	£12.30
BOT3228	2003-0004	125	50 x 101	24	12	£14.49
BOT3230	2003-0008	250	61 x 133	24	12	£20.90
BOT3232	2003-0016	500	73 x 170	28	12	£30.70
BOT3234	2003-0032	1000	91 x 216	38	6	£25.20