

Guardtack Mats

Cleanroom Entry Mats. Multi-layered, clean-film, full adhesive layers form mats that effectively capture dirt and dust from footborne traffic and equipment wheels before they enter your controlled environment. Each of the numbered layers of the mat peels off easily to reveal a new clean surface, 30 sheets per mat. The mats have a pressure sensitive adhesive that contains a spore growth retardant to ensure there can be no bacterial growth on the mat. The adhesive also ensures the mat is held firmly in place. Company logo and messages can be added to your mats - please ask for details.

- Peel-off sheets eliminate the need for messy cleaning and make it easy to maintain a clean surface
- Smooth, effective surface is not damaged by shoes, booties or cart wheels
- Numbered corner tabs for easy tracking of remaining sheets
- Composition: LDPE (Low Density Polyethylene)
- Thickness/sheet: 0.05mm \pm 5% thickness/mat: 1.50mm \pm 5%
- Peel adhesion: 0.15kN/m; medium tack
- Colours: Blue, white, grey
- 30 Layers per mat



Code	Dims, inches	Colour	Pack
CLE9038	18 x 22 inches	White	4
CLE9003	18 x 45 inches	Blue	4
CLE9005	18 x 45 inches	Grey	4
CLE9007	18 x 45 inches	White	8
CLE9009	18 x 45 inches	Blue	8
CLE9012	18 x 45 inches	Grey	8
CLE9001	18 x 45 inches	White	4
CLE9046	24 x 36 inches	Blue	4
CLE8990	26 x 45 inches	White	4
CLE8992	26 x 45 inches	Blue	4
CLE8994	26 x 45 inches	Grey	4
CLE8996	26 x 45 inches	White	8
CLE8998	26 x 45 inches	Blue	8
CLE9014	36 x 36 inches	White	4
CLE9016	36 x 36 inches	Blue	4
CLE9020	36 x 36 inches	White	8
CLE9022	36 x 36 inches	Blue	8
CLE9024	36 x 36 inches	Grey	8
CLE9032	36 x 45 inches	White	4
CLE9034	36 x 45 inches	Blue	4

Clean Room Mats

Disposable mats for use without frame. Colour white. Low ion content. Simply peel off back strip and lay mat where required. Low profile, only 2mm thick. Size 1143 x 495mm. 30 layers per mat.

Code	Description	Pack
CLE4030	Disposable mats	4

