

Analox O₂NE+ Oxygen Depletion Safety Monitor Continued

Specifications	
Range, % O ₂	0.1 to 25
IP Rating	IP65 (IP43 alarm repeater)
Sounder, dBA	85 @ 10cm
Monitor Display	4 digit LCD in %
Warm Up Time	30 seconds
Operating Temperature, °C	0 to 40
Monitor Dims, mm	175 x 110 x 75
Alarm Repeater Dims, mm	72 x 170 x 45
Accuracy (5 to 25%), % O ₂	Better than ±0.75

Code	Alt Ref	Description	Price
GAS1302	AX1BK10X55Q1Y01	O ₂ NE+ oxygen depletion safety monitor	£882.00

Accessories

Code	Alt Ref	Description	Price
GAS1304	REP1QXX01	Additional alarm repeater for entrances	£94.00
GAS1306	SA70SREGKIT	Regulator kit	£332.00
GAS1308	9100-2650	Spare O ₂ sensor	£451.00
GAS1322	SA70S+CAL	Calibration kit	£575.00
GAS1324	8000-0069A	Flow adapter	£63.00



5

YEARS

5 year warranty on vacuum

Statebourne Cryolab LN₂ Storage Dewars



Aluminium liquid storage dewars, specifically designed for storing and decanting liquid nitrogen in a range of capacities from 2 to 50L. They combine the qualities of high thermal efficiency, easy handling and maximum safety. Ultra low evaporation rates. These are complemented by a range of accessories to allow for safety when handling liquid gases.

Code	CRY1000	CRY1002	CRY1004	CRY1006	CRY1008	CRY1010	CRY1012
LN ₂ Capacity, L	2	5	10	20	25	35	50
Evaporation Rate, L/day	0.08	0.18	0.2	0.18	0.25	0.27	0.28
Static Hold Time, days	25	28	50	111	100	130	179
Weight (empty), kg	2.8	4.1	5.6	9.1	7.4	12.2	15.5
Weight (full), kg	4.4	8.1	13.7	25.3	27.6	40.5	55.9
Neck Ø, mm	32	51	51	51	51	63	63
Overall Dims, h x Ø, mm	385 x 185	495 x 216	560 x 305	570 x 407	585 x 395	650 x 407	810 x 407

Code	Alt Ref	Model	Price
CRY1000	99-01015	Cryolab 2	£529.00
CRY1002	99-01030	Cryolab 5	£576.00
CRY1004	99-01040	Cryolab 10	£608.00
CRY1006	99-01050	Cryolab 20	£825.00
CRY1008	9901062	Cryolab 25	£785.00
CRY1010	99-01080	Cryolab 35	£873.00
CRY1012	99-01090	Cryolab 50	£1087.00

Accessories

Dispensing packages available, please enquire for more details.

Code	Alt Ref	Description	For use with	Price
CRY1100	9901064	Tipping trolley, 1033 x 472mm (h x w)	All Cryolab	£402.00
CRY1101	97-1057	Withdrawal device	Cryolab 35, 50	£689.00
CRY1102	99-01063	Roller base, 110mm (h)	Cryolab 25	£306.00
CRY1104	99-02082	Roller base, 110mm (h)	Cryolab 20, 35, 50	£306.00