


Gel Materials



Reagents and materials for electrophoretic gels.

Code	Description	Pack	Price
17046801	Agarose IEF	10g	£165.00
17132901	PlusOne Bromophenol blue	10g	£50.00
17131401	CHAPS	1g	£62.00
17600318	DeStreak reagent	1mL	£48.70
17600319	DeStreak rehydration solution	3mL	£105.00
17131801	Dithiothreitol (DTT)	1g	£39.70
17131802	Dithiothreitol (DTT)	5g	£113.00
17133501	Immobiline Drystrip cover fluid	1000mL	£68.00
RPN6302	Iodoacetamide	25g	£88.00
17130402	PlusOne N,N'-Methylene-bis-acrylamide	100g	£131.00
17130202	PlusOne Acrylamide PAGE	1000g	£134.00
17131101	PlusOne Ammonium persulphate	25g	£31.20
17133001	PlusOne Bind-Silane	25mL	£50.00
17132501	PlusOne Glycerol 85%	1000mL	£61.00
17132301	PlusOne Glycine	500g	£48.30
17131001	PlusOne ReadySol IEF 40% T 3% C	1000mL	£224.00
17133201	PlusOne Repel-Silane ES	500mL	£100.00
17131301	SDS (Sodium dodecyl sulphate)	100g	£56.00
17131201	PlusOne TEMED	25mL	£39.90
RPN6301	Thiourea	100g	£74.00
17132101	TRIS	500g	£60.00
17131901	Urea	500g	£59.00



Immobiline™ DryStrip Gels



Immobiline™ DryStrip Gels (IPG strips) are used for isoelectric focusing (IEF), run as the first dimension of 2-D electrophoresis or as a separate application, and are available in fourteen overlapping pH ranges and five strip lengths (7, 11, 13, 18, and 24cm).

- Maximise resolution, loading capacity, and reproducibility of first-dimension IEF
- Choose from a comprehensive range of overlapping gels that cover narrow (1pH unit), medium (3 - 5pH units), and broad (8pH units) pH ranges
- Select the broad pH 3 - 11 NL gradient gel for fast and efficient screening to gain a broad overview of total protein distribution
- Higher resolution can be obtained using narrow-range (1pH unit) gels
- Use matching IPG buffers for optimised running conditions
- Use IPG strips directly after reswelling without further manipulation
- Barcoded strips for convenient sample and batch tracking

Supplied in packs of 12

Code	Description	Dims, mm	Price	Code	Description	Dims, mm	Price
17600245	pH 3-10 NL	240	£137.00	17600374	pH 3-11 NL	110	£128.00
17600244	pH 3-10	240	£137.00	17600376	pH 3-11 NL	180	£117.00
17600377	pH 3-11 NL	240	£139.00	17600373	pH 3-11 NL	70	£126.00
17600246	pH 4-7	240	£139.00	18101660	pH 4-7	110	£95.00
17600247	pH 6-9	240	£372.00	17600113	pH 4-7	130	£110.00
17600115	pH 3-10 NL	130	£110.00	17123301	pH 4-7	180	£114.00
17123501	pH 3-10 NL	180	£112.00	17600110	pH 4-7	70	£77.00
17600112	pH 3-10 NL	70	£77.00	17600196	pH 6-11	130	£316.00
18101661	pH 3-10	110	£95.00	17600197	pH 6-11	180	£170.00
17600114	pH 3-10	130	£110.00	17600195	pH 6-11	110	£338.00
17123401	pH 3-10 L	180	£112.00	17600194	pH 6-11	70	£169.00
17600111	pH 3-10	70	£77.00	17600375	pH 3-11 NL	130	£155.00

