

## Whatman PTFE Membranes



Hydrophobic, chemically inert and suitable for applications involving aggressive organic solvents, strong acids, and alkalis. Applications include preparing samples for HPLC analysis, air and gas sterilisation in fermentation vessels and filtrate tanks, and the clarification of aggressive liquids.

Supplied in packs of 100

Code	Alt Ref	Ø, mm	Pore Size, µm	Price
FIL3102	7582-002	25	0.2	£334.00
FIL3112	7590-002	25	1	£339.00
FIL3104	7582-004	47	0.2	£429.00
FIL3108	7585-004	47	0.5	£439.00
FIL3114	7590-004	47	1	£460.00




## Whatman Nylon Filter Membranes



High-quality nylon membranes are suitable for filtering aqueous solutions, most organic solvents and are hydrophilic so eliminating the need for wetting agents. The membranes are flexible, durable and tear resistant, and can be autoclaved at 135°C. Suitable for use with a wide range of biological preparations, nylon membranes are often used where other membranes are unsuitable. White/plain.

Code	Alt Ref	Ø, mm	Pore Size, µm	Pack	Price
FIL3120	7402001	13	0.2	100	£101.00
FIL3128	7404001	13	0.45	100	£101.00
FIL3130	7404002	25	0.45	100	£114.00
FIL3124	7402004	47	0.2	100	£155.00
FIL3132	7404004	47	0.45	100	£155.00
FIL3138	7408004	47	0.8	100	£142.00
FIL3126	7402009	90	0.2	50	£175.00
FIL3134	7404009	90	0.45	50	£176.00



## Whatman Mixed Cellulose Ester (WME) Membranes



Whatman mixed cellulose ester membranes are 47mm in diameter, and offer a variety of gridded formats for microbial colony quantification. Composed of cellulose acetate and cellulose nitrate, they offer a high degree of internal surface area for greater product adsorption.

Supplied in packs of 100

Code	Alt Ref	Description	Sterility	Pore Size, µm	Price
FIL3028	7153-104	Black, gridded, with pad	Sterile	0.45	£209.00
FIL3058	7187-114	White, gridded, without pad	Sterile	0.2	£92.00
FIL3070	7141-004	White, black grid	Non sterile	0.45	£150.00
FIL3072	7140-104	White, plain, with pad	Non sterile	0.45	£159.00
FIL3074	7141-104	White, gridded, with pad	Sterile	0.45	£90.00
FIL3076	7141-114	White, gridded, without pad	Sterile	0.45	£90.00
FIL3078	7141-204	White, black grid, autoclavable	Non sterile	0.45	£101.00

