

## Epredia Shandon CTM6 Coverslipper

epredia | Authorized Distributor

Get compact automation of routine glass coverslipping. With its compact footprint the CTM6 is the ideal instrument to free staff from the time-consuming task of routine coverslipping.

- Operator safety - loading and unloading of racks without opening the hood
- Internal vapour extraction system with fan and permanent air-flow avoids any leakage of vapours at openings
- Integrated charcoal filtration
- Bubble-free mountant refill, 100mL reservoir
- Tool-free cleaning and maintenance
- Easy conversion to other coverslip sizes - all required parts included



Code	Alt Ref	Dims, w x d x h, mm	Price
HIS3602	970010	275 x 410 x 350	£25497.00

## Epredia Shandon ClearVue™ Coverslipper

epredia | Authorized Distributor

Boost speed and accuracy in busy laboratories by automating the coverslipping workload. Simply place a basket of slides onto the load rail and close the door - the ClearVue™ Coverslipper does all the rest.

- Accepts slides directly from the Thermo Scientific™ Varistain™ Gemini ES and Varistain™ 24-4, Sakura™ DRS2000 and Leica™ Auto-Stainer
- Permits loading of slide baskets from any of these stainers in a mix-or-match fashion without requiring operator adjustments to the coverslipper
- Automatically recognizes different types of slide preparations and delivers the correct amount of mountant for coverslipping each
- Handles both histology and cytology samples simultaneously without user interaction while still delivering right amount of mountant
- Handles up to 11 slide baskets at the same time with load-on-demand capabilities, accepts baskets loaded one at a time or up to five baskets loaded together
- Recognizes the position of each slide in a basket, removes each slide with a set of slide-grippers and returns every slide back, while maintaining slide positioning
- Intuitive user interface features a touchscreen panel
- Software monitors the type of slide preparation, size and number of slips dispensed and the number of baskets waiting to be coverslipped
- Supplied with user guide, 2 x ClearVue™ coverslip hoppers, 24 x 40mm, #1.5–1000 slips, 3 x ClearVue™ slide baskets (Gemini-compatible) with black sliders, 2 x mountant purge trays, ventilation filter, 2 x soft hair brushes and 3 x suction cups



Code	Alt Ref	Dims, l x w x h, mm	Weight, kg	Price
HIS3604	A79200001	575 x 645 x 500	50	£30344.00

### Accessories

Code	Alt Ref	Description	Pack	Price
HIS3606	A79210050	ClearVue™ slip hopper no. 1.5, 24 x 40mm	2 x 500	£81.00
HIS3608	A79210051	ClearVue™ slip hopper no. 1.5, 24 x 50mm	2 x 500	£69.00
HIS3610	A79210052	ClearVue™ slip hopper no. 1.5, 24 x 55mm	2 x 500	£79.00
HIS3612	A79210053	ClearVue™ slip hopper no. 1, 24 x 40mm	2 x 500	£91.00
HIS3614	A79210166	ClearVue™ slip hopper no. 1, 24 x 50mm	2 x 500	£92.00
HIS3616	A79210064	Gemini compatible baskets, black sliders	5	£102.00
HIS3618	A79210065	Gemini compatible baskets, white sliders	5	£108.00