



## SLS Lab Basics Variable Volume Pipettes

### Single Channel

- Ergonomic design ensuring light weight & soft plunger movement
- Click-stop system for volume selection
- Colour coded for easy identification
- Calibration confirms to EN-ISO 8655 with report enclosed with every pipette


Code	Channel	Volume Range, $\mu\text{L}$	Compatible Tips (SLS Lab Pro)	Price
SLS4900	Single	0.1 - 2.5	10 $\mu\text{L}$ /10XL	£115.00
SLS4902	Single	0.5 - 10	10 $\mu\text{L}$ /10XL	£115.00
SLS4910	Single	2 - 20	10XL	£115.00
SLS4904	Single	5 - 50	20L/100 $\mu\text{L}$	£115.00
SLS4906	Single	10 - 100	100 $\mu\text{L}$ /200 $\mu\text{L}$	£115.00
SLS4912	Single	20 - 200	200 $\mu\text{L}$ /300 $\mu\text{L}$	£115.00
SLS4908	Single	100 - 1000	1000 $\mu\text{L}$ /1250 $\mu\text{L}$	£115.00
SLS4914	Single	500 - 5000	Requires tips PIP7736	£138.00
SLS4916	Single	1000 - 10000	Requires tips PIP3640	£138.00



### Multichannel

- Smooth and light tip loading and ejection
- Leak free sealing of tips
- Manifold can rotate 360° for easy right or left hand operation
- Easily removable & fully autoclavable manifold

Code	Channel	Volume Range, $\mu\text{L}$	Compatible Tips (SLS Lab Pro)	Price
SLS4926	8	0.5 - 10	10 $\mu\text{L}$ /10XL	£520.00
SLS4918	8	5 - 50	20 $\mu\text{L}$ /100XL	£520.00
SLS4920	8	10 - 100	100 $\mu\text{L}$ /200 $\mu\text{L}$	£520.00
SLS4922	8	20 - 200	200 $\mu\text{L}$ /300 $\mu\text{L}$	£520.00
SLS4924	8	30 - 300	200 $\mu\text{L}$ /300 $\mu\text{L}$	£520.00
SLS4936	12	0.5 - 10	10 $\mu\text{L}$ /10XL	£367.00
SLS4928	12	5 - 50	20 $\mu\text{L}$ /100XL	£367.00
SLS4930	12	10 - 100	100 $\mu\text{L}$ /200 $\mu\text{L}$	£367.00
SLS4932	12	20 - 200	200 $\mu\text{L}$ /300 $\mu\text{L}$	£367.00
SLS4934	12	30 - 300	200 $\mu\text{L}$ /300 $\mu\text{L}$	£367.00



## KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



At a glance



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty



Sustainable product