



Sartorius Picus® Electronic Pipettes

SARTORIUS

Picus offers all the advantages you would expect from a Sartorius pipette whilst the NxT version is targeted at customers working in a regulated environment where full compliance is a necessity.

Picus® Electronic Pipettes

Incorporating an electronic brake and piston control system help ensure accurate and precise pipetting results whilst an intuitive user interface with an adjustment wheel ensure fast volume setting and easy menu navigation.

- Electronic tip ejection
- Built in plate tracker functionality guides pipetting steps for improved microplate work
- Light at 100g
- Optiload ensures tip loading with a constant and equal force for perfect sealing
- Hot key access for up to 10 saved programmes

Code	Alt Ref	Channel	Volume Range, μL	Price
PIP3800	735021	Single	0.2 - 10	£543.00
PIP3802	735041	Single	5 - 120	£543.00
PIP3804	735061	Single	10 - 300	£543.00
PIP3806	735081	Single	50 - 1000	£543.00
PIP3808	735101	Single	100 - 5000	£543.00
PIP3810	735111	Single	500 - 10000	£543.00
PIP3812	735321	8	0.2 - 10	£838.00
PIP3814	735341	8	5 - 120	£838.00
PIP3816	735361	8	10 - 300	£838.00
PIP3818	735391	8	50 - 1200	£838.00
PIP3820	735421	12	0.2 - 10	£1003.00
PIP3822	735441	12	5 - 120	£1003.00
PIP3824	735461	12	10 - 300	£1003.00
PIP3826	735491	12	50 - 1200	£1003.00

Picus® NxT Electronic Pipettes

Offers all the advantages of the Picus® electronic pipettor with the added benefits of RFID capability, password protection and repeat blowout. Models covering the volume range 0.2-10000 μL .

- Single, 8 and 12 channel versions available all with electronic tip ejection

Code	Alt Ref	Channel	Volume Range, μL	Price
PIP3900	LH-745021	Single	0.2 - 10	£639.00
PIP3901	LH-745041	Single	5 - 120	£639.00
PIP3902	LH-745061	Single	10 - 300	£639.00
PIP3903	LH-745081	Single	50 - 1000	£639.00
PIP3904	LH-745101	Single	100 - 5000	£639.00
PIP3905	LH-745111	Single	500 - 10000	£639.00
PIP3906	LH-745321	8	0.2 - 10	£986.00
PIP3907	LH-745341	8	5 - 120	£986.00
PIP3908	LH-745361	8	10 - 300	£986.00
PIP3909	LH-745391	8	50 - 1200	£986.00
PIP3910	LH-745421	12	0.2 - 10	£1169.00
PIP3911	LH-745441	12	5 - 120	£1169.00
PIP3912	LH-745461	12	10 - 300	£1169.00
PIP3913	LH-745491	12	50 - 1200	£1169.00

Clarity in Servicing

Our service options:

- Fault Find & Fix**
- Service Contracts**
- On-Site Repairs**
- Pipette Servicing**



To book your service or for further information please contact us...

Email:
servicecentre@scientific-labs.com

Telephone:
0115 977 5075

Website:
www.scientificlabs.co.uk/servicecentre

**Our Knowledge + Your Equipment
= Service Tailored to You**