

Suretech Biogel® Latex Gloves

Unsupported natural rubber latex glove, sterile gamma irradiated and powder free, with micro-rough finish. Offering an excellent dexterity, easy donning, and exceptional levels of comfort due to the inner Biogel® technology. Packaged in individually sealed film pair packs.

- Length approx 290mm, thickness (double wall) approx 0.50mm
- In compliance with the standards EN420 (General requirements), EN388 (Mechanical hazards) and EN374 (Chemical and micro-organisms hazards) by the notified laboratory BSI
- Manufactured under EC quality assurance system carried out by the BSI

Supplied in packs of 50

Code	Size	Price
SAF4021	6.5	£162.00
SAF4022	7	£162.00
SAF4023	7.5	£162.00
SAF4024	8	£162.00
SAF4025	8.5	£162.00



GlovePlus Vinyl Examination Gloves



GlovePlus powder-free vinyl is ideal for use in areas of low risk where the purpose of wearing a glove is to keep the user's hands clean. Soft, durable and powder free.

- EN455 Parts 1 & 2, AQL 1.5
- Free from all REACH listed chemicals, 100% NR Latex free
- Below 2mg/g glove powder content, DOP [DEHP] free
- Length: 240mm
- Duration of use: 0-30 minutes
- Not suitable for contact with bloodborne pathogens or chemicals
- Risk category low

Supplied in packs of 100

Code	Alt Ref	Size	Price
SAF5648	SV32	Small	£12.07
SAF5650	SV33	Medium	£12.07
SAF5652	SV34	Large	£12.07
SAF5654	SV35	X large	£12.07



uvex Profapren CF33 Chemical Protection Gloves



High quality unsupported polychloroprene black chemical safety glove for protection against a broad spectrum of chemicals and solvents. Chloroprene glove with flocked cotton liner.

- Good grip in damp and wet areas thanks to the grip structure in the palm
- Silicone free
- Good combination of flexibility and strength
- Tactile feel
- EN399 (3131) EN ISO 374-1:2016/Type A (AKLMNO)
- Length: 330mm

Supplied in packs of 10 pairs

Code	Alt Ref	Description	Size	Price
SAF8388	60119-07	Small	7	£34.30
SAF8390	60119-08	Medium	8	£34.30
SAF8392	60119-09	Large	9	£34.30
SAF8394	60119-10	X large	10	£34.30

