



Jenway 73 Series Spectrophotometers

Jenway

A range of single beam spectrophotometers which use icon software and have an improved navigation system. The 7305 and 7315 models use a press to read xenon lamp to provide accurate readings and extend lamplife. The 7310 and 7315 models include measurements modes for photometrics, concentration, spectrum scanning, quantitation and kinetics.

- All models have photometrics and concentration modes
- Bandwidth: 5nm
- Icon driven software
- Results and method saving to USB memory stick
- Autologging capabilities
- Press to read xenon lamp (7315)

Code	Alt Ref	Model	Description	Wavelength Range, nm	Light Source	Price
SPE4600	730001	7300	Visible, single beam	320 - 1000	Tungsten halogen	£1966.00
SPE4620	730501	7305	Visible/UV, single beam	198 - 1000	Xenon	£2949.00
SPE4640	731001	7310	Visible, scanning, single beam	320 - 1000	Tungsten halogen	£2676.00
SPE4660	731501	7315	Visible/UV, scanning, single beam	198 - 1000	Xenon	£3560.00

Accessories

Jenway

Code	Alt Ref	Description	Price
SPE4668	660101	Internal printer	£433.00
SPE4670	735401	Automatic 8 position cuvette holder	£1101.00
SPE4672	735201	Sipper pump	£1207.00
SPE4674	735301	Peltier	£823.00
SPE4676	735701	Combined sipper/peltier	£1615.00
SPE3514	630204	10 x 10mm path length cuvette holder	£136.00
SPE3518	630005	Adjustable path length cuvette holder 10 to 100mm	£136.00
SPR4682	637071	16/24mm test tube holder	£220.00
SPE3520	630304	Micro cuvette holder with adjustable aperture	£136.00
SPR4686	736201	Water heated 10 x 10mm path length cuvette holder	£489.00
SPR4688	735601	Boiling tube holder	£220.00
SPR4690	735001	Dust cover	£80.00

Take a closer look at our microscope range!

Including monocular, binocular
and trinocular models.

See pages **628-655** or visit
www.scientificlabs.co.uk

