

Stuart Overhead Stirrers Continued




SS20 High Performance Overhead Stirrer

A powerful stirrer designed for demanding laboratory applications. Suitable for stirring liquids up to medium viscosity, including oils and microbiological media, up to 25L capacity. Features a hollow shaft arrangement to facilitate easy adjustment of paddle height. Advanced overload protection and sophisticated, responsive speed control.

- Used with volumes up to 25L
- Low to medium viscosity liquids
- Hollow shaft
- Keyless chuck


Code	Alt Ref	Description	Price
STI2333	SS20	Overhead stirrer, high performance	£880.00



Stands

Extra strong retort stands with H-pattern base for stability and robust support rod. Both versions include bosshead. Heavy duty stand recommended for use with SS30 overhead stirrer.

Code	Alt Ref	Description	Rod Dims, Ø x L, mm	Price
STI2334	SS10/1	Stand, general purpose	16 x 700	£264.00
STI2337	SS10/2	Stand, heavy duty	25 x 850	£688.00



Stirrer Paddles

Made from high grade stainless steel, there is a choice of five paddle heads and two lengths of 8mm diameter paddle rod. The heads screw on and off the rods, so they can be mixed and matched depending upon stirring requirements. This gives maximum versatility and value for money.

Code	Alt Ref	Description	Price
STI2340	SS10/5	Paddle rod, 350mm	£23.10
STI2342	SS10/6	Paddle rod, 550mm	£28.40
STI2344	SS10/10	Small paddle head	£36.80
STI2346	SS10/11	Large paddle head	£34.70
STI2348	SS10/12	Large propeller head (2 blade)	£36.80
STI2350	SS10/13	Small propeller head (4 blade)	£34.70
STI2352	SS10/14	Anchor paddle head	£34.70