



TitroLine® TL7500 KF / TL7000

Titration Offers for SLS Customers

AUTOMATIC TITRATION MAKING MOISTURE CONTENT ANALYSIS EASY

The TitroLine® 7500 KF for volumetric KF titration

The TL7500 from SI Analytics makes moisture content analysis easy. With one press of the touch-sensitive buttons the HD screen gives you step-by-step instructions on how to do your analysis.

With the capability to attach batch numbers, store the info to USB stick and even output to 21 CFR part 11 software, the TL7500 KF is the ideal instrument for small industrial / food labs right up to major pharmaceuticals.

FREE Installation worth up to £1,000

Purchase any TitroLine® TL7500 KF titrator between May and December 2021 and receive FREE installation worth in the region of £500 - £1,000 (reagents supplied separately).

Standard terms and conditions apply. Offer only valid to mainland UK & Ireland only. Offer ends December 30th 2021. Subject to availability. See website for details.



Official Distributor

SI Analytics
a xylem brand

Titration Offers for SLS Customers AVAILABLE THROUGHOUT 2021



TitroLine® TL7000 Titrator

Containing over 50 built-in methods the TitroLine® TL7000 titrator is a great choice for every lab wishing to automate titrations whilst also having the ability to input new methods at any point in the future.

With the flexibility of an exchangeable head and standard electrode connection sockets, changing between titration types/methods is fast and easy with the TL7000.

The need to re-validate new methods is negated by allowing users to input their own mathematics if and when required.

The TL7000 is the ideal choice for every laboratory, with the capability to attach autosamplers and compliance with 21 CFR part 11.

Claim a FREE set of pH buffers

Purchase any TL7000 package and receive 1 set of pH buffers

Standard terms and conditions apply. Offer only valid to mainland UK & Ireland only. Offer ends December 30th 2021. Subject to availability. See website for details.

TitroLine® Titrators from SI Analytics, a Xylem brand



TitroLine® TL7500 KF Titrator

- Polarisable electrode input for set endpoint titrations („Dead-stop“ + KF)
- With standard methods for various applications such as titer determination, sample titration with 1- or 2-component reagent, etc.
- Volumetric KF titrations
- Also suitable for dosing tasks and sample preparation
- Formula editor for individual calculations
- User management/password



TitroLine® TL7000 Titrator

- Reagent data is securely stored in the intelligent and interchangeable modules
- High resolution pH/mV-electrode and temperature inputs for pH, ISE, redox (ORP) or photometric titrations.
- Perfect for non-aqueous titration such as acid and base numbers in oils (TAN and TBN)
- Export the results as PDF or CSV, also to networks.

xylem
Let's Solve Water

© 2021 Xylem. All rights reserved. SI Analytics is a trademark of Xylem or one of its subsidiaries.

Official Distributor

SI Analytics
a xylem brand

☎ Lisburn: 028 9262 1100
☎ Dublin: 01 4069464

@ sales@analab.co.uk or sales@analab.ie
www www.analab.co.uk or www.analab.ie

