

Motic AE2000 Trionocular Microscope Package

The AE2000 inverted microscope is the ideal instrument for routine live cell inspection. Superb optical quality can be expected from a brand new optical series of both brightfield and phase contrast lenses, following Motic's CCIS® Infinity optical system.

- ✓ Eyepieces: N-WF10X/20mm with diopter adjustment, +/-5 diopter
- ✓ Objectives: 4X/0.10 (WD 12.6mm), PH10X/0.25 (WD 4.1mm), LWD PH20X/0.30 (WD 4.7mm), LWD 40X/0.50 (WD 3.0mm)
- ✓ Package includes: AE2000 Trinocular Microscope, 0.5X C-mount camera adapter, Motacam 1080 BMH with integrated 1080p HD LCD screen (no computer required)



Worth over €1875

Free Camera
with
Mount and
Integrated
Screen

Code	Description	Offer Price, £	Offer Price, €
F9196	Motic AE2000 Trionocular Microscope Package	3103.00	3599.00